



piccolo^{xpress}[®]

chemistry analyzer

FAST, EASY, ACCURATE ON-SITE CHEMISTRY TESTING

REAL-TIME BLOOD CHEMISTRY DIAGNOSTICS FOR URGENT CARE CENTERS

BE PREPARED FOR WHATEVER WALKS IN YOUR DOOR

- On-site diagnostics for rapid patient triage
- Quick, accurate and informative patient management
- Treat the patients you can and confidently refer others with the appropriate diagnostic information
- Comprehensive care increases patient satisfaction and repeat visits
- Maximize profitability of each patient visit
- Provide greater service & value for the community you serve

DIAGNOSE AND TREAT A WIDE VARIETY OF CONDITIONS

- Abdominal pain – liver, renal, gall bladder and pancreatic function
- Metabolic disorders – vomiting, diarrhea and dizziness
- Weakness and fatigue – dehydration, electrolyte imbalances
- Hyperlipidemia and hypertension
- Medication initiation and monitoring

○ FAST AND DECISIVE

From draw to printed results in 3 easy steps and approximately 12 minutes



1 / Add sample



2 / Insert disc



3 / Read results

○ SIMPLE AND SMART

Fully automated system requires no special operating skills

- Self contained reagent discs deliver complete results in real time
- Very small sample – only 100 microliters of whole blood, serum or plasma
- As easy to use as a CD player – load sample, insert disc and read results
- Barcoded discs self calibrate with every run and minimize data input errors
- CLIA waived – any of your clinicians can use it
- Comprehensive chemistry results in a single use disc
- Portable so that you can take it wherever it's needed
- 5 bi-directional USB ports for easy connectivity to EMR/LIS or peripherals

○ ACCURATE AND DEPENDABLE

Amazing technology ensures accurate results every time

- Precision and specificity comparable to expensive high-throughput analyzers
- Long 12-month shelf life ensures minimal waste
- Use discs directly from the refrigerator – no warming required
- iQC[®] system provides continuous quality control checks on every run to eliminate the potential for erroneous results
- Panels are fully reimbursable through medicare and private payors



3240 Whipple Rd., Union City, CA 94587 | 800.822.2947 | www.piccoloxpress.com

©2015 Abaxis, Inc., Piccolo Xpress[®] and iQC[®] are registered trademarks of Abaxis, Inc. 888-3221 Rev.C