



i-STAT Alinity v Cartridge Test Menu

		CG4+	CG8+	G	Crea	6+	CHEM8+	EC8+
Hematology	Hematocrit (Hct)		•			•	•	•
	Hemoglobin (Hb)*		•			•	•	•
Chemistry	Blood Urea Nitrogen (BUN)					•	•	•
	Creatinine (Crea)				•		•	
	Ionized Calcium (iCa)		•				•	
	Glucose (Glu)		•	•		•	•	•
Electrolytes	Chloride (Cl)					•	•	•
	Sodium (Na)		•			•	•	•
	Potassium (K)		•			•	•	•
Acid Base	pH	•	•					•
	Partial Pressure of Carbon Dioxide (PCO ₂)	•	•					•
	Bicarbonate (HCO ₃)*	•	•					•
	Total Carbon Dioxide (TCO ₂)*	•	•				•	•
	Anion Gap (AnGap)*						•	•
	Base Excess (BE)*	•	•					•
Blood Gas	Partial Pressure of Oxygen (PO ₂)	•	•					
	Oxygen Saturation (sO ₂)*	•	•					
Specialty	Lactate (Lac)	•						

*Calculated Value

To place your order today, contact your local authorized Abaxis Distributor, or for additional information contact Abaxis Customer Service:

800.822.2947

www.abaxis.com

