






Analytes	ALB	ALP	ALT	AMY	AST	BA	BUN	Ca	CHOL	CHW	CK	Cl-	CRE	GGT	GLOB*	GLU	K+	Mg	Na+	PHB	PHOS	T4	TBIL	tCO <sub>2</sub>	TP	UA
<b>Canine Wellness Profile #500-0044</b> 	•	•	•				•	•		•			•		•	•					•		•		•	
<b>Comprehensive Diagnostic Profile #500-0038</b>	•	•	•	•			•	•					•		•	•	•		•		•		•		•	
<b>Preventive Care Profile Plus #500-0047</b> 	•	•	•		•		•	•				•	•		•	•	•		•				•	•	•	
<b>Prep Profile II #500-0026</b>		•	•				•						•			•									•	
<b>Electrolyte Plus #500-0046</b> 												•					•		•					•		
<b>Critical Care Plus #500-0042</b>			•				•						•	•		•	•		•				•			
<b>Kidney Profile Plus #500-0045</b> 	•						•	•				•	•			•	•		•		•		•		•	
<b>T4/Cholesterol Profile #500-0037</b>									•													•				
<b>Phenobarbital Profile #500-0049</b> 	•	•	•		•		•							•						•			•			
<b>Mammalian Liver Profile #500-0040</b>	•	•	•			•	•		•					•									•			
<b>Equine Profile Plus #500-0043</b>	•				•		•	•			•		•	•	•	•	•		•				•	•	•	
<b>Avian/Reptilian Profile Plus #500-0041</b>	•				•	•	•				•				•	•	•		•		•				•	•
<b>Large Animal Profile #500-0023</b>	•	•			•		•	•			•			•	•			•		•					•	

\*Calculated value



For VetScan VS2 Only

FOR VETERINARY USE ONLY

To place your order today, contact your local authorized Abaxis Distributor, or for additional information contact Abaxis Customer Service:

 800.822.2947

 [www.abaxis.com](http://www.abaxis.com)

 [vetscan@abaxis.com](mailto:vetscan@abaxis.com)