


# VETSCAN<sup>®</sup> Phenobarbital Profile



Phenobarbital is an anticonvulsant drug used to prevent seizures in dogs and cats. Monitoring the level of phenobarbital is essential to achieve maximal seizure control while avoiding harmful side effects. Overdosing can lead to liver damage so it is crucial to carefully dose phenobarbital and closely monitor phenobarbital levels to ensure a safe and effective treatment.

The VetScan Phenobarbital Profile is a cost-effective rotor that evaluates the level of phenobarbital, and provides additional liver values all in one panel.

Featured Tests:



- Phenobarbital (PHB)
- Alkaline phosphatase (ALP)
- Alanine aminotransferase (ALT)
- Aspartate aminotransferase (AST)
- Gamma Glutamyl Transferase (GGT)
- Total bilirubin (TBIL)
- Albumin (ALB)
- Blood Urea Nitrogen (BUN)



## VetScan VS2 Chemistry Analyzer:

The Phenobarbital Profile is available exclusively for use with the Abaxis VetScan VS2 Chemistry Analyzer. **It's just one more example of practice-changing innovation that comes standard with every VS2.**

- Small sample size - only 2 drops of whole blood, serum or plasma
- Speed: 12 minutes to results
- Ease of use
- Small footprint and portable
- Little to no maintenance
- Intelligent Quality Control (iQC) on every run

## Intelligent Quality Control (iQC)

With the VetScan VS2 Chemistry Analyzer, quality control is built-in—thanks to a sophisticated iQC system that combines with proprietary algorithms to assure quality and dependability with every result. No more worrying about running controls and calibrations.



# VETSCAN® Chemistry Comprehensive Test Portfolio

NEW TEST

Analytes	ALB	ALP	ALT	AMY	AST	BA	BUN	Ca	CHOL	CHW	CK	Cl-	CRE	GGT	GLOB*	GLU	K+	Mg	Na+	PHB	PHOS	T4	TBIL	tCO <sub>2</sub>	TP	UA
<b>Canine Wellness Profile #500-0044</b>	•	•	•				•	•		•			•		•	•					•		•		•	
<b>Comprehensive Diagnostic Profile #500-0038</b>	•	•	•	•			•	•					•		•	•	•		•		•		•		•	
<b>Preventive Care Profile Plus #500-0047</b>	•	•	•		•		•	•				•	•		•	•	•		•				•	•	•	
<b>Prep Profile II #500-0026</b>		•	•				•						•			•									•	
<b>Electrolyte Plus #500-0046</b>												•					•		•					•		
<b>Critical Care Plus #500-0042</b>			•				•					•	•			•	•		•					•		
<b>Kidney Profile Plus #500-0045</b>	•						•	•				•	•			•	•		•		•			•		
<b>T4/Cholesterol Profile #500-0037</b>									•													•				
<b>Phenobarbital Profile #500-0049</b>	•	•	•		•		•							•						•			•			
<b>Mammalian Liver Profile #500-0040</b>	•	•	•			•	•		•					•									•			
<b>Equine Profile Plus #500-0043</b>	•				•		•	•			•		•	•	•	•	•		•				•	•	•	
<b>Avian/Reptilian Profile Plus #500-0041</b>	•				•	•		•			•				•	•	•		•		•				•	•
<b>Large Animal Profile #500-0023</b>	•	•			•		•	•			•			•	•			•		•					•	

\*Calculated value



For VetScan VS2 Only

FOR VETERINARY USE ONLY

To place your order today, contact your local authorized Abaxis Distributor, or for additional information contact Abaxis Customer Service:

800.822.2947

www.abaxis.com

vetscan@abaxis.com