


VETSCAN® Preventive Care Profile Plus



Obtain preventive care testing to identify subclinical diseases for patients at any age with the VetScan Preventive Care Profile Plus.

The Preventive Care Profile Plus is an advanced panel for pre-surgical, general health, ill patient, geriatric, and wellness testing that includes full electrolyte values and early information regarding acid-base status.

Featured Tests:



With few hands-on steps and no extra training required, every user can run 15 chemistries critical to monitoring their patient's needs with accurate results in just 12 minutes. **The Preventive Care Profile Plus is the most complete rotor profile from Abaxis, with 15 analytes readily available in one profile.**



VetScan VS2 Chemistry Analyzer:

The Preventive Care Profile Plus is available exclusively for use with the Abaxis VetScan VS2 Chemistry Analyzer.

It's just one more example of practice-changing innovation that comes standard with every VS2.

- Small sample size - only 2 drops of whole blood, serum or plasma
- Speed: 12 minutes to results
- Ease of use
- Small footprint and portable
- Little to no maintenance
- Intelligent Quality Control (iQC) on every run

Intelligent Quality Control (iQC)

With the VetScan VS2 Chemistry Analyzer, quality control is built-in—thanks to a sophisticated iQC system that combines with proprietary algorithms to assure quality and dependability with every result. No more worrying about running controls and calibrations.



VETSCAN® Chemistry Comprehensive Test Portfolio

NEW TEST

Analytes	ALB	ALP	ALT	AMY	AST	BA	BUN	Ca	CHOL	CHW	CK	Cl-	CRE	GGT	GLOB*	GLU	K+	Mg	Na+	PHB	PHOS	T4	TBIL	tCO ₂	TP	UA
Canine Wellness Profile #500-0044	•	•	•				•	•		•			•		•	•					•		•		•	
Comprehensive Diagnostic Profile #500-0038	•	•	•	•			•	•					•		•	•	•		•		•		•		•	
Preventive Care Profile Plus #500-0047	•	•	•		•		•	•				•	•		•	•	•		•				•	•	•	
Prep Profile II #500-0026		•	•				•						•			•									•	
Electrolyte Plus #500-0046												•					•		•					•		
Critical Care Plus #500-0042			•				•					•	•			•	•		•					•		
Kidney Profile Plus #500-0045	•						•	•				•	•			•	•		•		•			•		
T4/Cholesterol Profile #500-0037									•													•				
Phenobarbital Profile #500-0049	•	•	•		•		•							•						•			•			
Mammalian Liver Profile #500-0040	•	•	•			•	•		•					•									•			
Equine Profile Plus #500-0043	•				•		•	•			•		•	•	•	•	•		•				•	•	•	
Avian/Reptilian Profile Plus #500-0041	•				•	•		•			•				•	•	•		•		•				•	•
Large Animal Profile #500-0023	•	•			•		•	•			•			•	•			•		•					•	

*Calculated value



For VetScan VS2 Only

FOR VETERINARY USE ONLY

To place your order today, contact your local authorized Abaxis Distributor, or for additional information contact Abaxis Customer Service:

800.822.2947

www.abaxis.com

vetscan@abaxis.com