



		CG4+	CG8+	G	Crea	E3+	6+	CHEM8+	EC8+	ACT Celite®	cTnl
Hematology	Hematocrit (Hct)		•			•	•	•	•		
	Hemoglobin (Hgb)		•			•	•	•	•		
Chemistry	Blood Urea Nitrogen (BUN)						•	•	•		
	Creatinine				•			•			
	Ionized Calcium (iCa)		•					•			
	Glucose (Glu)		•	•			•	•	•		
Electrolytes	Chloride (Cl)*					•	•	•	•		
	Sodium (Na)		•			•	•	•	•		
	Potassium (K)		•			•	•	•	•		
Acid Base	pH	•	•						•		
	PCO ₂	•	•						•		
	HCO ₃	•	•						•		
	TCO ₂	•	•					•	•		
	Anion Gap							•	•		
	Base Excess	•	•						•		
Blood Gas	PO ₂	•	•								
	sO ₂	•	•								
Specialty	Lactate	•									
	ACT Celite									•	
	Cardiac Troponin I (cTnl)										•

*Chloride test on the E3+ cartridge is only available on the VetScan i-STAT 1 analyzer.

To place your order today, contact your authorized Abaxis Distributor or for additional information contact Abaxis Customer Service at 1.800.822.2947. www.abaxis.com