

11 OF 13 AVAILABLE PANELS ARE CLIA WAIVED

The most comprehensive CLIA waived menu available.

**Comprehensive Metabolic Panel**

ALB, ALP, ALT, AST, BUN, Ca, Cl<sup>-</sup>, CRE, GLU, K<sup>+</sup>, Na<sup>+</sup>, TBIL, tCO<sub>2</sub>, TP

**Basic Metabolic Panel**

BUN, Ca, Cl<sup>-</sup>, CRE, GLU, K<sup>+</sup>, Na<sup>+</sup>, tCO<sub>2</sub>

**Lipid Panel**

CHOL, CHOL/HDL\*, HDL, LDL\*, TRIG, VLDL\* \*Calculated

**Lipid Panel Plus**

ALT, AST, CHOL, CHOL/HDL\*, GLU, HDL, LDL\*, TRIG, VLDL\*

**Liver Panel Plus**

ALB, ALP, ALT, AMY, AST, GGT, TBIL, TP

**General Chemistry 13**

ALB, ALP, ALT, AMY, AST, BUN, Ca, CRE, GGT, GLU, TBIL, TP, UA

**General Chemistry 6**

ALT, AST, BUN, CRE, GGT, GLU

**Electrolyte Panel**

Cl<sup>-</sup>, K<sup>+</sup>, Na<sup>+</sup>, tCO<sub>2</sub>

**Kidney Check**

BUN, CRE

**MetLyte 8 Panel**

BUN, CK, Cl<sup>-</sup>, CRE, GLU, K<sup>+</sup>, Na<sup>+</sup>, tCO<sub>2</sub>

**Renal Function Panel**

ALB, BUN, Ca, Cl<sup>-</sup>, CRE, GLU, K<sup>+</sup>, Na<sup>+</sup>, PHOS, tCO<sub>2</sub>



**Basic Metabolic Panel Plus**

BUN, Ca, Cl<sup>-</sup>, CRE, GLU, K<sup>+</sup>, Na<sup>+</sup>, tCO<sub>2</sub>, Mg, Lactate Dehydrogenase

**Hepatic Function Panel**

ALB, ALP, ALT, AST, DBIL, TBIL, TP



COMPREHENSIVE CLIA WAIVED CHEMISTRY  
FOR HEALTH SCREENING OPERATIONS

piccolo *xpress*  
chemistry analyzer

ABAXIS

# Health Screening Just Became Comprehensive!

**piccolo** *xpress*<sup>™</sup>  
chemistry analyzer

With the Piccolo Xpress<sup>™</sup> chemistry analyzer, you can now conduct comprehensive CLIA Waived screenings on-site in minutes.

Health screeners, wellness and occupational health organizations traditionally rely on limited screening capabilities at the point-of-care. Lipid and glucose screenings alone do little to provide the complete health picture. And, managing a multitude of samples to be screened at the lab, then following up with patients afterward is an ineffective, time-consuming and often expensive process.

With the Piccolo Xpress chemistry analyzer in your arsenal, patients can now have on-site laboratory results for critical functions such as:

- Liver
- Intestine
- Parathyroid
- Renal
- Pancreas
- Bone Marrow
- Metabolic
- Respiratory
- Jaundice
- Cardiovascular
- Hepatic
- Shock
- Gall bladder
- Diabetic
- and much more

## Three simple steps is all it takes.



## Taking health screening to a new level

Health screening operations play a vital role in promoting wellness and early detection of health problems. By adding the Piccolo Xpress chemistry analyzer, health screening operations can:

- Conduct comprehensive labs on-site
- Reduce costs while increasing value of offering
- Gain a competitive advantage over rival providers
- Maximize operational efficiency and effectiveness
- Increase patient and customer satisfaction from complete service

## Exactly what the doctor ordered

When a physician performs a well-person check-up, they order the same battery of tests that are available through the Piccolo Xpress. The Piccolo allows health screening clinicians to:

- Deploy a comprehensive CLIA waived menu of tests
- Receive lab-accurate results in 10-12 minutes
- Provide physicians with complete, timely and accurate health analysis in advance of appointments
- Employ advanced technology that's easy to use
- Utilize only 100 microliters of whole blood
- Help patients improve overall health and wellness

## Designed with Health Screeners in mind

The Piccolo Xpress chemistry analyzer is perfectly suited for health screening operations. Its small and compact size fits nicely in tight spaces. Simple 3-step operation makes it easy for screeners to obtain accurate results. And, it's fast - taking only 30 seconds to prepare a sample and 10-12 minutes for a printed result. Walk-away operation means the Piccolo is self-sufficient and frees the screener's time to perform other patient examinations.

## Exceptionally Easy

The Piccolo Xpress's fully automated system requires no special operating skills

- 0.1cc whole blood, serum or plasma
- Self calibrates with every test run
- Built-in iQC<sup>™</sup> eliminates daily control testing and meets CLIA and CLSI guidelines
- Full-color touch screen with intuitive user interface
- Results are printed on a self-adhesive hard copy or via electronic transfer, allowing immediate access to critical information at the point-of-care
- Large memory capacity stores up to 5000 patient records
- Externally accessible bi-directional input/output and USB data transfer technology
- CLIA waived reagent discs eliminate the need for moderately complex lab requirements

## Exceptionally Accurate

Advanced technology ensures sound results

- Centrifugal action mixes the reagents and sample in a disposable self-contained disc – reducing procedural errors
- Accuracy and precision comparable to reference grade laboratory analyzers
- iQC<sup>™</sup> system performs continuous quality control checks on every run
- Every self-contained reagent disc is barcoded for automatic self-calibration



Interested? Contact Abaxis at: 1-800-822-2947 [www.abaxis.com](http://www.abaxis.com)

**ABAXIS**