


i-STAT® Level 2 Control

LOT 111098

 **Exp.** 2019-02-28

CLEW: A35

CG4+ & G3+ LOT D		\bar{x} (Mean)	R (Range)
pH		7.489	7.439 - 7.539
PCO2	mmHg	28.9	21.9 - 35.9
	kPa	3.85	2.92 - 4.79
PO2	mmHg	118	100 - 135
	kPa	15.7	13.3 - 18.0
Lactate/Lac	mmol/L	2.18	1.68 - 2.68
	mg/dL	19.6	15.1 - 24.1
	g/L	0.196	0.151 - 0.241

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748727-01

20-Dec-2017

i-STAT® Level 2 Control

LOT 111098

 **Exp.** 2019-02-28

CLEW: A35

CG8+, EG7+, EG6+, G3+, Crea & CG4+		\bar{x} (Mean)	R (Range)
LOT	A, K, L, M, N, W, Y		
Na	mmol/L, mEq/L	136	131 - 140
K	mmol/L, mEq/L	3.7	3.4 - 4.0
iCa	mmol/L	1.29	1.21 - 1.38
	mg/dL	5.2	4.8 - 5.5
	mEq/L	2.6	2.4 - 2.8
pH		7.462	7.412 - 7.512
PCO2	mmHg	30.5	23.5 - 37.5
	kPa	4.07	3.13 - 5.00
PO2	mmHg	117	99 - 134
	kPa	15.6	13.2 - 17.9
Glucose/Glu	mg/dL	122	106 - 139
	g/L	1.22	1.06 - 1.39
	mmol/L	6.8	5.9 - 7.7
Creatinine/Crea	mg/dL	1.7	1.3 - 2.1
	umol/L	150	115 - 186
Lactate/Lac	mmol/L	2.16	1.66 - 2.66
	mg/dL	19.5	15.0 - 24.0
	g/L	0.195	0.150 - 0.240

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748727-01

20-Dec-2017

i-STAT® Level 2 Control

LOT 111098

 **Exp.** 2019-02-28

CLEW: A35

CHEM8+ LOT H		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	135	130 - 139
K	mmol/L, mEq/L	3.7	3.4 - 4.0
Cl	mmol/L, mEq/L	90	84 - 96
iCa	mmol/L	1.30	1.21 - 1.38
	mg/dL	5.2	4.8 - 5.5
	mEq/L	2.6	2.4 - 2.8
TCO2	mmol/L, mEq/L	24	18 - 31
Glucose/Glu	mg/dL	124	107 - 140
	g/L	1.24	1.07 - 1.40
	mmol/L	6.9	5.9 - 7.8
Creatinine/Crea	mg/dL	1.8	1.4 - 2.2
	umol/L	159	124 - 194
BUN	mg/dL	10	7 - 13
Urea	mmol/L	3.6	2.5 - 4.6
	mg/dL	21	15 - 28
	g/L	0.21	0.15 - 0.28

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748727-01

20-Dec-2017

i-STAT® Level 2 Control

LOT 111098

Exp. 2019-02-28

CLEW: A35

EC8+, 6+, EC4+, E3+ & G LOT F, H, J, K		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	136	131 - 140
K	mmol/L, mEq/L	3.7	3.4 - 4.0
Cl	mmol/L, mEq/L	89	83 - 95
pH		7.393	7.343 - 7.443
PCO2	mmHg	30.3	23.3 - 37.3
	kPa	4.04	3.11 - 4.97
Glucose/Glu	mg/dL	121	105 - 138
	g/L	1.21	1.05 - 1.38
	mmol/L	6.7	5.8 - 7.7
BUN	mg/dL	10	7 - 13
Urea	mmol/L	3.6	2.5 - 4.6
	mg/dL	21	15 - 28
	g/L	0.21	0.15 - 0.28

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748727-01

20-Dec-2017