

i-STAT® Level 2 Control

LOT 111094

 **Exp.** 2018-10-31

CLEW: A35

CG4+ & G3+ LOT D		\bar{x} (Mean)	R (Range)
pH		7.484	7.434 - 7.534
PCO2	mmHg	29.1	22.1 - 36.1
	kPa	3.88	2.95 - 4.81
PO2	mmHg	118	100 - 135
	kPa	15.7	13.3 - 18.0
Lactate/Lac	mmol/L	2.23	1.73 - 2.73
	mg/dL	20.1	15.6 - 24.6
	g/L	0.201	0.156 - 0.246

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748301-01

31-Oct-2017

i-STAT® Level 2 Control

LOT 111094

 **Exp.** 2018-10-31

CLEW: A35

CG8+, EG7+, EG6+, G3+, Crea & CG4+ LOT A, K, L, M, N, W, Y		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	135	131 - 140
K	mmol/L, mEq/L	3.7	3.4 - 4.0
iCa	mmol/L	1.30	1.21 - 1.38
	mg/dL	5.2	4.8 - 5.5
	mEq/L	2.6	2.4 - 2.8
pH		7.456	7.406 - 7.506
PCO2	mmHg	30.8	23.8 - 37.8
	kPa	4.11	3.17 - 5.04
PO2	mmHg	116	99 - 134
	kPa	15.5	13.2 - 17.9
Glucose/Glu	mg/dL	122	105 - 138
	g/L	1.22	1.05 - 1.38
	mmol/L	6.8	5.8 - 7.7
Creatinine/Crea	mg/dL	1.7	1.3 - 2.1
	umol/L	150	115 - 186
Lactate/Lac	mmol/L	2.21	1.71 - 2.71
	mg/dL	19.9	15.4 - 24.4
	g/L	0.199	0.154 - 0.244

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748301-01

31-Oct-2017

i-STAT® Level 2 Control

LOT 111094

 **Exp.** 2018-10-31

CLEW: A35

CHEM8+ LOT H		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	135	130 - 139
K	mmol/L, mEq/L	3.7	3.4 - 4.0
Cl	mmol/L, mEq/L	90	84 - 96
iCa	mmol/L	1.31	1.22 - 1.39
	mg/dL	5.2	4.9 - 5.6
	mEq/L	2.6	2.4 - 2.8
TCO2	mmol/L, mEq/L	24	18 - 31
Glucose/Glu	mg/dL	123	106 - 140
	g/L	1.23	1.06 - 1.40
	mmol/L	6.8	5.9 - 7.8
Creatinine/Crea	mg/dL	1.7	1.3 - 2.1
	umol/L	150	115 - 186
BUN	mg/dL	10	7 - 13
Urea	mmol/L	3.6	2.5 - 4.6
	mg/dL	21	15 - 28
	g/L	0.21	0.15 - 0.28

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748301-01

31-Oct-2017

i-STAT® Level 2 Control

LOT 111094

 **Exp.** 2018-10-31

CLEW: A35

EC8+, 6+, EC4+, E3+ & G		\bar{x} (Mean)	R (Range)
LOT	F, H, J, K		
Na	mmol/L, mEq/L	136	131 - 140
K	mmol/L, mEq/L	3.7	3.4 - 4.0
Cl	mmol/L, mEq/L	89	83 - 95
pH		7.388	7.338 - 7.438
PCO2	mmHg	30.7	23.7 - 37.7
	kPa	4.09	3.16 - 5.03
Glucose/Glu	mg/dL	121	104 - 137
	g/L	1.21	1.04 - 1.37
	mmol/L	6.7	5.8 - 7.6
BUN	mg/dL	10	7 - 13
Urea	mmol/L	3.6	2.5 - 4.6
	mg/dL	21	15 - 28
	g/L	0.21	0.15 - 0.28

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748301-01

31-Oct-2017