

FERRET REFERENCE RANGES FOR COMMON ANALYTES

For VetScan® and VetScan® VS2

ANALYTE	RANGE	UNITS	N
ALBUMIN	1.9 – 3.8	g/dL	98
ALP	8 - 72	IU/L	97
ALT	65 - 346	IU/L	98
AMYLASE	4 - 50	IU/L	96
BILE ACIDS (BA)	1 - 8	µmol/L	103
BUN (Urea Nitrogen)	10 - 38	mg/dL	100
CALCIUM	8.0 – 10.4	mg/dL	99
CHOLESTEROL	102 - 245	mg/dL	104
CREATININE	0.2 – 0.7	mg/dL	98
GGT	5 – 15	IU/L	103
GLU	65 - 145	mg/dL	99
GLOBULIN	2.3 – 4.5	g/dL	98
PHOSPHOROUS	3.6 – 7.3	mg/dL	98
POTASSIUM	4.1 – 5.5	mmol/L	99
SODIUM	146 - 156	mmol/L	98
TOTAL BILIRUBIN	0.3 – 0.6	mg/dL	100
TOTAL PROTEIN	5.0 – 7.6	g/dL	100

Courtesy of Karen Rosenthal, DVM, University of Pennsylvania Veterinary Hospital, Philadelphia, PA

REVISED 12/2009